



## Use case 1: Power On

1. **Name:** Power On
2. **Actor:** User
3. **Pre-conditions:** AED is turned off
4. **Main Success Scenario:**
  - a. User presses the “Power On” button
  - b. A pop-up window appears, informing the user that the machine has done a self-check, and everything is operational.
  - c. AED audio plays, confirming ready for use.
5. **Extensions:**
  - a. If the AED encounters issues during the self-check(step 4b), an error message is provided and guidance for resolution.
6. **Success Guaranteed:** AED is powered on and the user receives confirmation of successful complete self-check and operational status.




## Use case 2: Pad Already Placed on Patient

1. **Name:** Pad Placement
2. **Actor:** User
3. **Pre-conditions:** AED is turned on and passes self-check
4. **Main Success Scenario:**
  - a. User selects the “Electrode Placement” option
  - b. AED provides a simulated prompt, guiding the user through the process of placing the pads accurately on the patient's bare chest.
5. **Extensions:**
  - a. If the user places the pad incorrectly(step 4b), the AED provides guidance for correct placement.
  - b. If the user is still struggling with pad placements(step 4b), the AED offers a step-by-step instruction until the pads are placed properly.
  - c. If there is a loss of connection during the process, it goes back to placing pads.
6. **Success Guaranteed:** Electrode pads are placed correctly on the patients chest, to ensure accurate readings.



### **Use case 3: Regular Heartbeat Reading(PEA)**

1. **Name:** Regular Heartbeat Reading(PEA)
2. **Actor:** AED simulation
3. **Pre-conditions:** AED is turned on and electrode pads are placed properly
4. **Main Success Scenario:**
  - a. AED analyzes the patient's heart rhythm and detects a regular heartbeat(PEA).
  - b. The AED provides audio instructions to monitor the patient while waiting for help.
5. **Extensions:**
  - a. If the AED detects a shockable rhythm(step 4a), it moves to the appropriate scenarios(e.g. irregular heartbeat VF or VT).
  - b. If there is a loss of connection during the process, it goes back to placing pads.
6. **Success Guaranteed:** AED recognizes the regular heartbeat and instructs the user to monitor the patient.




## Use case 4: Regular heartbeat (Flatlined)

1. **Name:** Regular Heartbeat (Flatlined)
2. **Actor:** AED Simulation
3. **Pre-conditions:** AED is turned on and electrode pads are placed properly.
4. **Main Success Scenario:**
  - a. AED analyzes the patient's heart rhythm and detects a flatline.
  - b. The AED provides audio and visual prompts indicating the need for CPR.
  - c. The AED displays helpful tips, guiding the user on the proper chest compression depth(2 to 4 inches) and providing real-time feedback.
  - d. If compression depth is insufficient, the AED displays “Bad compression” and guides the user on improvement technique(s).
  - e. After CPR initiation, the AED continues to analyze for shocks, providing real-time feedback and instructions.

## 5. Extensions:

- If the AED detects a shockable rhythm during analysis(step 4a), it transitions to the appropriate scenario( e.g. irregular heartbeat VF or VT).
- If there is a loss of connection during the process, it goes back to placing pads.


## 6. Success Guaranteed: AED recognizes flatlined rhythm, guides the user through CPR and provides feedback.



## Use case 5: Irregular heartbeat (Ventricular Fibrillation(VF))

- Name:** Irregular Heartbeat (VF)
- Actor:** AED Simulation
- Pre-conditions:** AED is turned on and electrode pads are placed properly.
- Main Success Scenario:**
  - AED analyzes the patient's heart rhythm and detects VF.
  - The AED provides audio and visual prompts indicating the availability of shocks.
  - User initiates shock delivery, AED administers 200j shock(if administering more than 1 shock, then the AED should administer 300j shock from there on).

- 6. Success Guaranteed:** AED recognizes VF, delivers shocks, guides the user through transition to CPR if stabilization occurs.




1. **Name:** Irregular Heartbeat (VT)
2. **Actor:** AED Simulation
3. **Pre-conditions:** AED is turned on and electrode pads are placed properly.
4. **Main Success Scenario:**
  - a. AED analyzes the patient's heart rhythm and detects VT.
  - b. The AED provides audio and visual prompts indicating the availability of shocks.
  - c. User initiates shock delivery, AED administers 200j shock(if administering more than 1 shock, then the AED should administer 300j shock from there on).

- d. The AED prompts the user to perform CPR.
- e. The AED continues to analyze for additional shocks, delivering randomized 1-3 shocks if needed.
- f. If the patient stabilizes, the AED transitions to CPR automatically.

**5. Extensions:**

- a. If the AED detects a non-shockable rhythm during analysis(step 4a), it transitions to the appropriate scenario( e.g. regular heartbeat flatlined).
- b. If there is a loss of connection during the process, it goes back to placing pads.

**6. Success Guaranteed:** AED recognizes VT, delivers shocks, guides the user through transition to CPR if stabilization occurs.



## Use case 7: Battery Too Low


1. **Name:** Battery Too Low
2. **Actor:** AED Simulation
3. **Pre-conditions:** AED is turned on with a low battery.
4. **Main Success Scenario:**
  - a. AED detects a low battery.
  - b. A pop-up window notifies the user about low battery and the inability to provide shock to the patient.

c. AED instructs the user to replace the battery.

**5. Extensions:**

- a. If the low battery situation still persists despite battery replacement attempts(step 4c), the AED provides an error message.
- b. If there is a loss of connection during the process, it goes back to placing pads.

**6. Success Guaranteed:** AED notifies the user of low battery and instructs them to replace battery for proper functionality.




## Use Case 8: Loss of Connection


1. **Name:** Loss of Connection
2. **Actor:** AED Simulation
3. **Pre-conditions:** AED is turned on, and electrode pads are connected to the patient.
4. **Main Success Scenario:**



- a. AED detects a loss of connection between the electrode pads and the patient.
  - b. A pop-up window notifies the user about the loss of connection.
  - c. AED waits for the user to correctly place pads on the patient.
5. **Extensions:** None
6. **Success Guaranteed:** AED recognizes the loss of connection and notifies the user.







ID	Requirements	Related Use Case	Fulfilled By	Test	Description
1	Design GUI resembling AED Plus's display	Use Case1- Power on	AEDSimulation.c pp and Scenario.cpp	Interface usability test	Ensure GUI is intuitive and displays necessary info
2	Simulate Cardiac arrhythmia detection	Use case 5- irregular heartbeat (VF) Use case6- irregular heartbeat (VT)	AEDSimulation.c pp and Scenario.cpp	Functional test	Detect VF and Vt rhythms; provide visual/textual prompts
3	Provide real-time CPR feedback	Use Case 3- regular heartbeat Reading. Use case4- Regular Heart(flatlined)	AEDSimulation.c pp and Scenario.cpp	Feedback accuracy test	Measure and assess chest compressions

<b>4</b>	Develop visual prompts for emergency scenarios	Usecase4- Regular Heartbeat(flatlined) Usecase5- Irregular Heartbeat (VF) Usecase6 – irregular heartbeat (VT)	AEDSimulation.cpp, Scenario.cpp and CPRFeedback.cpp	Scenario simulation test	Guide user through rescue situation
<b>5</b>	Interactive system for AED operation simulation	Use case1-power on Usecase2- pad placement Usecase3- Regular heartbeat reading. Usecase4- regular heartbeat(flat lined) Usecase5- irregular heartbeat Use case 6- irregular heartbeat Usecase7- battery too low	AEDSimulation.cpp , Scenario.cpp, Battery.cpp and CPRFeedback.cpp	User interaction test	Allow users to perform actions like in real world AED
<b>6</b>	Create simulated scenarios for software testing	Usecase5- irregular heartbeat (VF) Usecase6-irregular	AEDSimulation.cpp , Scenario.cpp and CPRFeedback.cpp	Scenario test	Test with various patient conditions

		heartbeat (VT)			
<b>7</b>	Display simulated ECG waveform and device status	Usecase1- power on Usecase2- pad placement Usecase3- regular heartbeat reading. Usecase4- regular heartbeat(flat lined) Usecase5-irr equalr heartbeat (VF) Usecase6- irregular heartbeat (VT) Usecase7-ba ttery too low	AEDSimulation.cpp and Battery.cpp	Display functionality test	Accurately simulate ECG and device status
<b>8</b>	Instruction for shock delivery and CPR in emergencies	Usecase4-R egular heartbeat(flat lined) Usecase5-irr equalr Use case5- irregular heartbeat heartbeat (VF) Usecase6-irr equalr heartbeat (VT0	AEDSimulation.cp, Scenario.cpp and CPRFeedback.cpp	Compliance test	Ensuring instructions are clear and accurate

	Pad Placement Check	Uses cases: 2,3,4,5,6,8.	Scenario.cpp		
<b>9</b>	Compression Depth Analysis	Use cases: 2,3,4,5,6,8	Scenario.cpp, CPRFeedback.cpp	Compliance and Scenario Test	The device is able to give the user feedback on their compressions, feedback on both the rate and depth of the compressions.
<b>10</b>	Device performs a self test	Use cases: 1,7	AEDSimulation.cpp, Scenario.cpp, Battery.cpp	Compliance and Scenario Test	The device should perform a self-test check, and lets the user know whether the battery is too low or if the pads are functional

### **Explanation of our design decision:**

**The mainwindow has event-listeners on the “Turn on” Button, when it is pressed it sends a request to the AEDSimulation class but that depends on the use case number selected from the drop down button (1-8)**

**The scenario class has a list of functions to call depending on the use case selected**

**When the heartbeat analysis is being conducted, the scenario class then sends a request to the AED Simulation class to perform an analysis**

**When the CPR is performed, the scenario class then sends a request to the AED Simulation class which then sends a request to the CPRFeedback class**

**Design Patterns used:**

- 1) Singleton pattern: Each of our classes have only one instance and they all provide a global access to them**
- 2) Observer pattern: When the different scenarios changes, all the other classes are changed accordingly**
- 3) Strategy Pattern: There are presence of switch statements which holds the function patterns depending on what the use case number is**